Cunhas Reflexivas

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PT Dispositivo
ES Médico

Translucid Wedges

The Translucid Wedges are indicated for general use, in restorative procedures using transparent matrices, with the mechanical function of adapting the matrix and/or separate adjacent teeth, and the optical function of helping the composite polymerization in the proximal region. The Translucid Wedges do not need clipping and/or adjustments and guarantee the matrix fixing, providing restorations without sub gingival excess and with perfect contour. Available in small and medium sizes.

Use: Select the desired wedge. Insert the wedge into the interproximal space. Once the procedure is completed, the wedge must be removed.

Composition: Polyvinyl chloride.

Warnings, precautions and restrictions: Professional use only. Single use only. Its reuse can cause cross contamination. If there is proven allergy to any of its components, the product must not be used. To avoid an accidental swallowing of this product, it is necessary the application of a total isolation technique of the operator field. The product has a sharp point and must be handled with caution to avoid risk of causing harm to the patient. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Storage: To be stored in a cool, protected from sunlight.

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